

| Designation ENG | Packaging                                      | Packaging | Shelf life | Minimal<br>Remaining<br>shelf life<br>(days) | Storing T° | Storage after preparation                | Autoclave      | Supplements |
|-----------------|------------------------------------------------|-----------|------------|----------------------------------------------|------------|------------------------------------------|----------------|-------------|
| <b>A1010GC</b>  | Bacteriological Agar type A                    | 5 kg      | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1010HA</b>  | Bacteriological Agar type A                    | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1012GC</b>  | Bacteriological Agar type E                    | 5 kg      | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1012HA</b>  | Bacteriological Agar type E                    | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1101HA</b>  | Bacteriological malt extract                   | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1202GC</b>  | Yeast extract                                  | 5 kg      | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1202HA</b>  | Yeast extract                                  | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1401GC</b>  | Tryptone USP                                   | 5 kg      | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1401HA</b>  | Tryptone USP                                   | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1402GC</b>  | Pancreatic Digest of Casein Codex              | 5 kg      | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1402HA</b>  | Pancreatic Digest of Casein Codex              | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1404HA</b>  | Acid Hydrolysate of Casein                     | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1601HA</b>  | Papaic digest of soybean meal USP              | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1702HA</b>  | Pancreatic digest of meat type 2               | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>A1710HA</b>  | Bacteriological meat extract                   | 500 g     | 5 years    | 455                                          | 2-30°C     | -                                        | -              | -           |
| <b>BK001HA</b>  | TSN Agar                                       | 500 g     | 5 years    | 455                                          | 2-30°C     | 0 days                                   | 15 min 121°C   | -           |
| <b>BK002HA</b>  | Brilliant Green Bile Broth (BGBB)              | 500 g     | 5 years    | 455                                          | 2-30°C     | 90 days, shielded from light             | 15 min 121°C   | -           |
| <b>BK003HA</b>  | Nutrient Broth                                 | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months 2-25°C                          | 15 min 121°C   | -           |
| <b>BK007HA</b>  | Chloramphenicol glucose Agar                   | 500 g     | 3 years    | 455                                          | 2-30°C     | vial 6 months, plate 30 days<br>at 2-8°C | 15 min 121°C   | -           |
| <b>BK010HA</b>  | Laurylsulfate-Tryptose Broth                   | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months 15-25°C, 0 days at<br>2-8°C     | 15 min 121°C   | -           |
| <b>BK011GC</b>  | VRBG Agar                                      | 5 kg      | 3 years    | 455                                          | 2-30°C     | 0 days                                   | NO             | -           |
| <b>BK011HA</b>  | VRBG Agar                                      | 500 g     | 3 years    | 455                                          | 2-30°C     | 0 days                                   | NO             | -           |
| <b>BK012HA</b>  | M 17 Broth                                     | 500 g     | 5 years    | 455                                          | 2-20°C     | 6 months 2-8°C                           | 20 min à 115°C | -           |
| <b>BK013HA</b>  | Wort Agar - Dehydrated base medium             | 500 g     | 5 years    | 455                                          | 2-30°C     | vial 6 months, plate 30 days<br>at 2-8°C | 15 min 110°C   | glycerol    |
| <b>BK014HA</b>  | Tryptone-salt Broth (Maximum Recovery Diluent) | 500 g     | 3 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                       | 15 min 121°C   | -           |
| <b>BK015HA</b>  | Brain-Heart Broth                              | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                       | 15 min 121°C   | -           |
| <b>BK017HA</b>  | Thioglycollate medium with resazurin           | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months at 2-8 °C                       | 15 min 121°C   | -           |
| <b>BK018GC</b>  | Buffered peptone water (25.5 g/L)              | 5 kg      | 5 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                       | 15 min 121°C   | -           |
| <b>BK018HA</b>  | Buffered peptone water (25.5 g/L)              | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                       | 15 min 121°C   | -           |

| Designation ENG | Packaging                              | Packaging | Shelf life | Minimal Remaining shelf life (days) | Storing T°                  | Storage after preparation                                        | Autoclave    | Supplements   |
|-----------------|----------------------------------------|-----------|------------|-------------------------------------|-----------------------------|------------------------------------------------------------------|--------------|---------------|
| <b>BK019HA</b>  | Columbia Agar - Dehydrated base medium | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-8°C                            | 15 min 121°C | sheep's blood |
| <b>BK020HA</b>  | CLED Agar                              | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-8°C                            | 20 min 115°C | -             |
| <b>BK022HA</b>  | SS Agar                                | 500 g     | 5 years    | 455                                 | 2-30°C                      | 15 days at 2-8°C shielded from light                             | NO           | -             |
| <b>BK023HA</b>  | Bromcresol Purple Lactose Agar         | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-8°C                            | 15 min 121°C | -             |
| <b>BK024HA</b>  | Meat liver glucose 0.6% Agar           | 500 g     | 5 years    | 455                                 | 2-30°C                      | 6 months at 2-8°C                                                | 15 min 121°C | -             |
| <b>BK025HA</b>  | Sabouraud Dextrose Agar                | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-25°C                           | 15 min 121°C | -             |
| <b>BK026HA</b>  | Sabouraud Broth                        | 500 g     | 4 years    | 455                                 | 2-30°C                      | 6 months at 2-25°C                                               | 15 min 121°C | -             |
| <b>BK027HA</b>  | Sabouraud Chloramphenicol Agar         | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-8°C                            | 15 min 121°C | -             |
| <b>BK028HA</b>  | Tryptone-Soy Agar (Blood Agar base)    | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-8°C                            | 15 min 121°C | sheep's blood |
| <b>BK029HA</b>  | Brain-Heart Agar                       | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-8°C                            | 15 min 121°C | -             |
| <b>BK030HA</b>  | Mannitol salt Agar                     | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-8°C                            | 15 min 121°C | -             |
| <b>BK031HA</b>  | TSC Agar - Dehydrated base medium      | 500 g     | 5 years    | 455                                 | 2-30°C                      | 6 months at 2-25°C                                               | 15 min 121°C | BS006         |
| <b>BK033HA</b>  | ROGOSA Agar                            | 500 g     | 5 years    | 455                                 | 2-20°C                      | not recommended                                                  | NO           | -             |
| <b>BK034HA</b>  | KLIGLER Agar                           | 500 g     | 5 years    | 455                                 | 2-30°C                      | non inclined tubes 6 months, tubes inclined 8 days at 2-8°C      | 15 min 121°C | -             |
| <b>BK036HA</b>  | DRIGALSKI lactose Agar                 | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial 6 months, plate 30 days at 2-8°C                            | 20 min 115°C | -             |
| <b>BK037HA</b>  | SLANETZ and BARTLEY Agar               | 500 g     | 5 years    | 455                                 | 2-30°C, shielded from light | vial not recommended, plate 30 days at 2-8°C shielded from light | NO           | -             |
| <b>BK040HA</b>  | TCBS Agar                              | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial not recommended, plate 8 days at 2-8°C                      | NO           | -             |
| <b>BK042HA</b>  | BCP Dextrose Agar                      | 500 g     | 4 years    | 455                                 | 2-30°C                      | 6 months at 2-8°C                                                | 15 min 121°C | -             |
| <b>BK045HA</b>  | Malt Agar                              | 500 g     | 5 years    | 455                                 | 2-30°C                      | vial/tube 6 months, plate 30 days at 2-8°C                       | 15 min 121°C | -             |

| Designation ENG | Packaging                                                  | Packaging | Shelf life | Minimal Remaining shelf life (days) | Storing T° | Storage after preparation                                       | Autoclave    | Supplements       |
|-----------------|------------------------------------------------------------|-----------|------------|-------------------------------------|------------|-----------------------------------------------------------------|--------------|-------------------|
| <b>BK046GC</b>  | Trypto-casein Soy Broth (TSB)                              | 5 kg      | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                              | 15 min 121°C | -                 |
| <b>BK046HA</b>  | Trypto-casein Soy Broth (TSB)                              | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                              | 15 min 121°C | -                 |
| <b>BK047GC</b>  | Trypto-casein Soy Agar (TSA)                               | 5 kg      | 5 years    | 455                                 | 2-30°C     | vial 6 months 2-25°C , plate 30 days at 2-8°C                   | 15 min 121°C | -                 |
| <b>BK047HA</b>  | Trypto-casein Soy Agar (TSA)                               | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months 2-25°C , plate 30 days at 2-8°C                   | 15 min 121°C | -                 |
| <b>BK048HA</b>  | MUELLER-HINTON Agar                                        | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months, plate 30 days at 2-8°C                           | 15 min 115°C | -                 |
| <b>BK049HA</b>  | Cetrimide Agar - Dehydrated base medium (without glycerol) | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months, plate 30 days at 2-8°C                           | 15 min 121°C | glycerol          |
| <b>BK050HA</b>  | MacCONKEY Agar                                             | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months, plate 30 days at 2-8°C                           | 15 min 121°C | -                 |
| <b>BK053HA</b>  | Oxytetracycline glucose Agar (OGA)                         | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months 2-25°C , plate 15 days at 2-8°C                   | 15 min 121°C | BS008             |
| <b>BK055GC</b>  | BAIRD-PARKER Agar - Dehydrated base medium                 | 5 kg      | 5 years    | 455                                 | 2-30°C     | base 6 months at 2-8°C, plate 30 days at 2-8°C                  | 15 min 121°C | 01BPRPF : BS045 c |
| <b>BK055HA</b>  | BAIRD-PARKER Agar - Dehydrated base medium                 | 500 g     | 5 years    | 455                                 | 2-30°C     | base 6 months at 2-8°C, plate 30 days at 2-8°C                  | 15 min 121°C | 01BPRPF : BS045 c |
| <b>BK056HA</b>  | EMB Agar (LEVINE)                                          | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months, plate 30 days at 2-8°C                           | 15 min 121°C | -                 |
| <b>BK057HA</b>  | ENDO Agar                                                  | 500 g     | 5 years    | 455                                 | 2-30°C     | vial not recommended, plate 3 days at 2-8°C shielded from light | 15 min 121°C | -                 |
| <b>BK058HA</b>  | XLD Agar                                                   | 500 g     | 3 years    | 455                                 | 2-30°C     | 8 days at 2-8°C                                                 | NO           | -                 |
| <b>BK060HA</b>  | Azide Dextrose Broth (acc. to ROTHE)                       | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-8°C                                               | 15 min 121°C | -                 |
| <b>BK061HA</b>  | Ethyl violet azide (LITSKY) Broth                          | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-8°C                                               | 15 min 121°C | -                 |
| <b>BK062HA</b>  | Desoxycholate (0.1%) lactose Agar                          | 500 g     | 5 years    | 455                                 | 2-30°C     | not recommended                                                 | NO           | -                 |
| <b>BK065HA</b>  | Desoxycholate lactose Agar                                 | 500 g     | 5 years    | 455                                 | 2-30°C     | not recommended                                                 | NO           | -                 |
| <b>BK067HA</b>  | Hektoen Enteric Agar                                       | 500 g     | 5 years    | 455                                 | 2-30°C     | vial not recommended, plate 30 days at 2-8°C                    | NO           | -                 |
| <b>BK070GC</b>  | MRS Broth                                                  | 5 kg      | 5 years    | 455                                 | 2-20°C     | 6 months at 2-8°C                                               | 15 min 121°C | -                 |
| <b>BK070HA</b>  | MRS Broth                                                  | 500 g     | 5 years    | 455                                 | 2-20°C     | 6 months at 2-8°C                                               | 15 min 121°C | -                 |
| <b>BK071HA</b>  | Brilliant Green Agar (KRISTENSEN)                          | 500 g     | 5 years    | 455                                 | 2-30°C     | vial not recommended, plate 30 days at 2-8°C                    | 15 min 121°C | -                 |
| <b>BK081HA</b>  | BUTTAUX-BROGNIART hypersalted Broth (double strength)      | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-8°C                                               | 15 min 121°C | -                 |

| Designation ENG | Packaging                                                      | Packaging | Shelf life | Minimal Remaining shelf life (days) | Storing T° | Storage after preparation                                 | Autoclave    | Supplements    |
|-----------------|----------------------------------------------------------------|-----------|------------|-------------------------------------|------------|-----------------------------------------------------------|--------------|----------------|
| <b>BK082HA</b>  | Lactose Broth                                                  | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                        | 15 min 121°C | -              |
| <b>BK084HA</b>  | Peptone water                                                  | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                        | 15 min 121°C | -              |
| <b>BK087HA</b>  | MAYEUX SANDINE ELLIKER Agar (MSE)                              | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months, plate 15 days at 2-8°C                     | 20 min 110°C | -              |
| <b>BK088HA</b>  | M 17 Agar                                                      | 500 g     | 5 years    | 455                                 | 2-20°C     | 6 months at 2-8°C                                         | 20 min 115°C | -              |
| <b>BK089HA</b>  | MRS Agar                                                       | 500 g     | 3 years    | 455                                 | 2-20°C     | 6 months at 2-8°C                                         | 15 min 121°C | -              |
| <b>BK090HA</b>  | Reinforced Clostridial Agar                                    | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-8°C                                         | 15 min 121°C | -              |
| <b>BK091HA</b>  | Brilliant Green Agar (EDEL and KAMPELMACHER)                   | 500 g     | 5 years    | 455                                 | 2-30°C     | 30 days at 2-8°C                                          | NO           | -              |
| <b>BK094HA</b>  | Reinforced Clostridial Medium (RCM) Broth                      | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-8°C                                         | 15 min 121°C | -              |
| <b>BK095HA</b>  | Potato Dextrose Agar (PDA)                                     | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months at 2-25°C, plate 30 days at 2-8°C           | 15 min 121°C | -              |
| <b>BK101HA</b>  | WILKINS-CHALGREN Agar                                          | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months, plate 30 days 2-8°C                        | 15 min 121°C | -              |
| <b>BK103HA</b>  | Orange Serum Agar (OSA)                                        | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months 2-25°C, plate 15 days at 2-8°C              | 15 min 115°C | -              |
| <b>BK107HA</b>  | MacCONKEY Broth purple                                         | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-8°C shielded from light                     | 15 min 121°C | -              |
| <b>BK108HA</b>  | MUELLER-HINTON Broth                                           | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                        | 15 min 121°C | -              |
| <b>BK112HA</b>  | <i>Listeria</i> enrichment Broth (LEB acc. to FDA)             | 500 g     | 5 years    | 455                                 | 2-30°C     | 30 days at 2-8°C shielded from light                      | 15 min 115°C | -              |
| <b>BK113HA</b>  | <i>Listeria</i> enrichment Broth (UVM I modified)              | 500 g     | 5 years    | 455                                 | 2-30°C     | 30 days at 2-8°C shielded from light                      | 15 min 115°C | -              |
| <b>BK115HA</b>  | FRASER Broth - Dehydrated base medium                          | 500 g     | 3 years    | 455                                 | 2-30°C     | Base 6 months, plate 30 days at 2-8°C shielded from light | 15 min 121°C | BS059 ou BS062 |
| <b>BK116HA</b>  | Bacillus cereus Agar (acc. to Mossel) - Dehydrated base medium | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months, plate 30 days 2-8°C                        | 15 min 121°C | BS055          |
| <b>BK118HA</b>  | CFC Agar - Dehydrated base medium                              | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months, plate 8 days 2-8°C                         | 15 min 121°C | BS022          |
| <b>BK119HA</b>  | Bromcresol Purple Lactose Broth                                | 500 g     | 5 years    | 455                                 | 2-30°C     | 30 days at 2-8°C                                          | 15 min 121°C | -              |
| <b>BK123HA</b>  | Tergitol 7 Agar - Dehydrated base medium                       | 500 g     | 3 years    | 455                                 | 2-30°C     | Base 6 months, plate 30 days at 2-8°C shielded from light | 15 min 121°C | BS026          |
| <b>BK126HA</b>  | Sugar-Free Agar (SFA)                                          | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months at 2-25°C, plate 30 days at 2-8°C           | 15 min 121°C | -              |

| Designation ENG | Packaging                                                         | Packaging | Shelf life | Minimal<br>Remaining<br>shelf life<br>(days) | Storing T° | Storage after preparation                                                                         | Autoclave     | Supplements         |
|-----------------|-------------------------------------------------------------------|-----------|------------|----------------------------------------------|------------|---------------------------------------------------------------------------------------------------|---------------|---------------------|
| <b>BK127HA</b>  | EE Broth (acc. to MOSSEL)                                         | 500 g     | 5 years    | 455                                          | 2-30°C     | 3 months at 2-8 °C shielded from light                                                            | 15 min 121°C  | -                   |
| <b>BK128HA</b>  | Buffered peptone Broth with sodium chloride pH 7.0                | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                                                                                | 15 min 121°C  | -                   |
| <b>BK129HA</b>  | SLANETZ and BARTLEY Agar - Dehydrated base medium                 | 500 g     | 5 years    | 455                                          | 2-30°C     | conservation si autoclavage vial 6 months, plate 30 days at 2-8°C                                 | NO (optional) | BS027               |
| <b>BK130HA</b>  | Count-tact Agar                                                   | 500 g     | 5 years    | 455                                          | 2-20°C     | vial 6 months at 2-25°C, plate 30 days at 2-8°C                                                   | 15 min 121°C  | -                   |
| <b>BK131GC</b>  | Buffered peptone water (20 g/L)                                   | 5 kg      | 5 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                                                                                | 15 min 121°C  | -                   |
| <b>BK131HA</b>  | Buffered peptone water (20 g/L)                                   | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                                                                                | 15 min 121°C  | -                   |
| <b>BK132HA</b>  | KF Streptococcus Agar - Dehydrated base medium                    | 500 g     | 5 years    | 455                                          | 2-30°C     | plate 8 days at 2-8°C                                                                             | NO            | BS027               |
| <b>BK133GC</b>  | FRASER Broth (II) - Dehydrated base medium                        | 5 kg      | 5 years    | 455                                          | 2-30°C     | FI base 6 months at 2-8°C, FI complete 6 months at 2-8°C (2 months at 2-25°C) shielded from light | 15 min 121°C  | BS030               |
| <b>BK133HA</b>  | FRASER Broth (II) - Dehydrated base medium                        | 500 g     | 5 years    | 455                                          | 2-30°C     | FI base 6 months at 2-8°C, FI complete 6 months at 2-8°C (2 months at 2-25°C) shielded from light | 15 min 121°C  | BS030               |
| <b>BK134HA</b>  | MSRV Agar - Dehydrated base medium                                | 500 g     | 5 years    | 455                                          | 2-30°C     | 8 days at 2-8°C shielded from light                                                               | NO            | BS033               |
| <b>BK135HA</b>  | MÜLLER-KAUFFMANN Broth - Dehydrated base medium                   | 500 g     | 5 years    | 455                                          | 2-30°C     | base 30 days 2-8°C shielded from light, complete not recommended                                  | NO            | tassium iodide + io |
| <b>BK137HA</b>  | Eugon LT 100 Broth                                                | 500 g     | 5 years    | 455                                          | 2-20°C     | 6 months at 2-25°C                                                                                | 15 min 121°C  | -                   |
| <b>BK138HA</b>  | Eugon LT 100 Agar                                                 | 500 g     | 5 years    | 455                                          | 2-20°C     | vial 6 months 2-25°C, plate 30 days 2-8°C                                                         | 15 min 121°C  | -                   |
| <b>BK140HA</b>  | Lactose-sulfite (LS) Broth                                        | 500 g     | 3 years    | 455                                          | 2-30°C     | not recommended                                                                                   | 15 min 121°C  | -                   |
| <b>BK141GC</b>  | BRYANT and BURKEY Broth with resazurin (modified acc. to BERGERE) | 5 kg      | 5 years    | 455                                          | 2-30°C     | not recommended                                                                                   | 15 min 121°C  | -                   |
| <b>BK141HA</b>  | BRYANT and BURKEY Broth with resazurin (modified acc. to BERGERE) | 500 g     | 5 years    | 455                                          | 2-30°C     | not recommended                                                                                   | 15 min 121°C  | -                   |
| <b>BK144GC</b>  | Plate count Agar (PCA)                                            | 5 kg      | 5 years    | 455                                          | 2-30°C     | vial 6 months 2-25°C, plate 30 days at 2-8°C                                                      | 15 min 121°C  | -                   |

| Designation ENG | Packaging                                             | Packaging | Shelf life | Minimal<br>Remaining<br>shelf life<br>(days) | Storing T° | Storage after preparation                                 | Autoclave    | Supplements         |
|-----------------|-------------------------------------------------------|-----------|------------|----------------------------------------------|------------|-----------------------------------------------------------|--------------|---------------------|
| <b>BK144HA</b>  | Plate count Agar (PCA)                                | 500 g     | 5 years    | 455                                          | 2-30°C     | vial 6 months 2-25°C, plate 30 days at 2-8°C              | 15 min 121°C | -                   |
| <b>BK145HA</b>  | PALCAM Agar - Dehydrated base medium                  | 500 g     | 5 years    | 455                                          | 2-30°C     | vial base 6 months, plate 30 days at 2-8°C                | 15 min 121°C | BS004               |
| <b>BK146HA</b>  | TBX Agar                                              | 500 g     | 5 years    | 455                                          | 2-20°C     | vial 6 months, plate 15 days at 2-8°C shielded from light | 15 min 121°C | -                   |
| <b>BK146HM</b>  | TBX Agar                                              | 100 g     | 5 years    | 455                                          | 2-20°C     | vial 6 months, plate 15 days at 2-8°C shielded from light | 15 min 121°C | -                   |
| <b>BK147HA</b>  | MacCONKEY-sorbitol (CT-SMAC) Agar - Dehydrated base   | 500 g     | 5 years    | 455                                          | 2-30°C     | vial base 6 months, plate 30 days at 2-8°C                | 15 min 121°C | BS037               |
| <b>BK148HA</b>  | RAPPAPORT-VASSILIADIS Soja (RVS) Broth                | 500 g     | 5 years    | 455                                          | 2-30°C     | tube 6 months at 2-8°C                                    | 15 min 115°C | -                   |
| <b>BK150HA</b>  | Modified Tryptone soy Broth (mTSB) - Dehydrated base  | 500 g     | 5 years    | 455                                          | 2-30°C     | base 6 months, complete 30 days at 2-8°C                  | 15 min 121°C | BS033               |
| <b>BK151HA</b>  | Rose Bengal Chloramphenicol Agar                      | 500 g     | 5 years    | 455                                          | 2-30°C     | vial 3 months, plate 7days at 2-8°C                       | 15 min 121°C | -                   |
| <b>BK152GC</b>  | VRBL Agar                                             | 5 kg      | 5 years    | 455                                          | 2-30°C     | not recommended                                           | NO           | -                   |
| <b>BK152HA</b>  | VRBL Agar                                             | 500 g     | 5 years    | 455                                          | 2-30°C     | not recommended                                           | NO           | -                   |
| <b>BK153HA</b>  | Yeast extract Agar                                    | 500 g     | 5 years    | 455                                          | 2-30°C     | fl 6 months 2-25°C, plate 30 days at 2-8°C                | 15 min 121°C | -                   |
| <b>BK157HA</b>  | Meat-liver glucose Agar                               | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months at 2-8°C                                         | 15 min 121°C | -                   |
| <b>BK158HA</b>  | BEA Agar                                              | 500 g     | 5 years    | 455                                          | 2-30°C     | fl 6 months, plate 30 days at 2-8°C                       | 15 min 121°C | -                   |
| <b>BK159HA</b>  | GIOLITTI and CANTONI Broth with Tween®80 - Dehydrated | 500 g     | 5 years    | 455                                          | 2-30°C     | base 6 months 2-8°C, complete not recommended             | 15 min 121°C | Potassium tellurite |
| <b>BK161HA</b>  | Plate count Agar with skimmed milk                    | 500 g     | 5 years    | 455                                          | 2-30°C     | vial 6 months at 2-25°C, plate 30 days at 2-8°C           | 15 min 121°C | -                   |
| <b>BK162HA</b>  | EC Broth                                              | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                                        | 15 min 121°C | -                   |
| <b>BK163HA</b>  | Tryptophan Broth                                      | 500 g     | 5 years    | 455                                          | 2-30°C     | 6 months at 2-25°C                                        | 15 min 121°C | -                   |
| <b>BK165HA</b>  | CN Agar for Pseudomonas - Dehydrated base medium (K)  | 500g      | 5 years    | 455                                          | 2-30°C     | Base 3 months, complete 30 days at 2-8°C                  | 15 min 121°C | glycerol            |
| <b>BK168HA</b>  | XLD Agar (ISO 6579-1)                                 | 500 g     | 5 years    | 455                                          | 2-30°C     | 8 days at 2-8°C                                           | NO           | -                   |



| Designation ENG | Packaging                                              | Packaging | Shelf life | Minimal Remaining shelf life (days) | Storing T° | Storage after preparation                                         | Autoclave                 | Supplements        |
|-----------------|--------------------------------------------------------|-----------|------------|-------------------------------------|------------|-------------------------------------------------------------------|---------------------------|--------------------|
| <b>BK169HA</b>  | MKTTn Broth - Dehydrated base medium                   | 500 g     | 5 years    | 455                                 | 2-30°C     | base 30 days, complete 8 days 2-8°C                               | NO                        | ium iodide + iode+ |
| <b>BK170HA</b>  | Dichloran-Glycerol (DG18) Agar - Dehydrated base medi  | 500 g     | 5 years    | 455                                 | 2-20°C     | vial 3 months, plate 8days at 2-8°C                               | 15 min 121°C              | glycerol           |
| <b>BK173GC</b>  | Half-FRASER Broth - Dehydrated base medium             | 5 kg      | 3 years    | 455                                 | 2-30°C     | 6 months 2-8°C, (or complete 2 months 2-25°C) shielded from light | 15 min 121°C              | BS059 ou BS062     |
| <b>BK173HA</b>  | Half-FRASER Broth - Dehydrated base medium             | 500 g     | 3 years    | 455                                 | 2-30°C     | 6 months 2-8°C, (orcomplete 2 months 2-25°C) shielded from light  | 15 min 121°C              | BS059 ou BS062     |
| <b>BK176HA</b>  | MRS Broth (Non-animal)                                 | 500 g     | 5 years    | 455                                 | 2-20°C     | 6 months 2-8°C                                                    | 15 min 121°C              | -                  |
| <b>BK179HA</b>  | R2A                                                    | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months 2-25°C, plate 30 days at 2-8°C                      | 15 min 121°C              | -                  |
| <b>BK183HA</b>  | COMPASS® Enterococcus Agar                             | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months plate 30 days at 2-8°C shielded from light          | 15 min 121°C              | -                  |
| <b>BK185HA</b>  | Nutrient Agar 2%                                       | 500 g     | 5 years    | 455                                 | 2-30°C     | vial 6 months at 2-25°C, plate 30 days at 2-8 °C                  | 15 min 121°C              | -                  |
| <b>BK186HA</b>  | Glutamate Broth                                        | 500 g     | 3 years    | 455                                 | 2-20°C     | 3 months at 2-8°C shielded from light                             | 15 min 115°C 10 min 116°C | -                  |
| <b>BK189HA</b>  | COMPASS® Bacillus cereus Agar                          | 500 g     | 3 years    | 455                                 | 2-30°C     | vial 6 months, plate 30 days at 2-8°C                             | 15 min 121°C              | supplement         |
| <b>BK191HA</b>  | MSRV Agar (ISO 6579-1)                                 | 500 g     | 3 years    | 455                                 | 2-30°C     | 15 days at 2-8°C shielded from light                              | NO                        | -                  |
| <b>BK192GC</b>  | COMPASS® Listeria Agar - Dehydrated base medium        | 5 kg      | 3 years    | 455                                 | 2-30°C     | vial 6 months (complete 4 h 44-47°C), plate 15 days at 2-8°C      | 15 min 121°C              | BS067 + BS068      |
| <b>BK192HA</b>  | COMPASS® Listeria Agar - Dehydrated base medium        | 500 g     | 3 years    | 455                                 | 2-30°C     | vial 6 months (complete 4 h 44-47°C), plate 15 days at 2-8°C      | 15 min 121°C              | BS067 + BS068      |
| <b>BK194GC</b>  | Salmonella Enrichment                                  | 5 kg      | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                                | 15 min 121°C              | -                  |
| <b>BK194HA</b>  | Salmonella Enrichment                                  | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                                | 15 min 121°C              | -                  |
| <b>BK195HA</b>  | SESAME Salmonella Detection                            | 500 g     | 3 years    | 455                                 | 2-30°C     | vial 6 months, plate 15 days at 2-8°C                             | NO                        | -                  |
| <b>BK198HA</b>  | Dichloran Rose Bengal Chloramphenicol (DRBC) Agar      | 500 g     | 3 years    | 455                                 | 2-8°C      | vial 3 months, plate 30 days at 2-8°C                             | 15 min 121°C              | -                  |
| <b>BK200HA</b>  | Cronobacter Chromogenic Isolation (CCI) agar - Dehydra | 500 g     | 3 years    | 455                                 | 2-30°C     | -                                                                 | -                         | -                  |

| Designation ENG | Packaging                                                  | Packaging | Shelf life | Minimal Remaining shelf life (days) | Storing T° | Storage after preparation                                                                  | Autoclave    | Supplements          |
|-----------------|------------------------------------------------------------|-----------|------------|-------------------------------------|------------|--------------------------------------------------------------------------------------------|--------------|----------------------|
| <b>BK202HA</b>  | COMPASS® Ecc Agar - Dehydrated medium                      | 500 g     | 3 years    | 455                                 | 2-30°C     | 3 months 2-8°C shielded from light                                                         | 15 min 121°C | -                    |
| <b>BK204HA</b>  | Chromogenic Coliform Agar (CCA)                            | 500 g     | 3 years    | 455                                 | 2-30°C     | 3 months 2-8°C shielded from light                                                         | 15 min 121°C | -                    |
| <b>BK205GC</b>  | Non-animal malt extract broth                              | 5 kg      | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                                                         | 15 min 121°C | -                    |
| <b>BK205HA</b>  | Non-animal malt extract broth                              | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                                                         | 15 min 121°C | -                    |
| <b>BK206GC</b>  | Thioglycollate medium with resazurin - Ready-to-use medium | 5 kg      | 5 years    | 455                                 | 2-30°C     | -                                                                                          | 15 min 121°C | -                    |
| <b>BK208HA</b>  | MKTTn Broth - Dehydrated base medium with novobiocin       | 500 g     | 3 years    | 455                                 | 2-30°C     | base 30 days, complete 8 days at 2-8°C                                                     | NO           | iodide and biosublin |
| <b>BK210HA</b>  | COMPASS® cc Agar                                           | 500 g     | 3 years    | 455                                 | 2-30°C     | base 6 months, complete plate 30 days 2-8°C shielded from light                            | 15 min 121°C | BS084                |
| <b>BK212HA</b>  | IRIS Salmonella Agar                                       | 500 g     | 4 years    | 455                                 | 2-8°C      | plate 30 days 2-8°C shielded from light                                                    | NO           | -                    |
| <b>BK214HA</b>  | RHAPSODY®- Dehydrated base medium                          | 500 g     | 4 years    | 455                                 | 2-30°C     | 15 days at 2-8°C shielded from light                                                       | -            | -                    |
| <b>BK215GC</b>  | Complete Half Fraser broth                                 | 5 kg      | 5 years    | 455                                 | 2-30°C     | -                                                                                          | -            | -                    |
| <b>BK215HA</b>  | Complete Half fraser broth                                 | 500 g     | 5 years    | 455                                 | 2-30°C     | -                                                                                          | -            | -                    |
| <b>BK216GC</b>  | EASY Staph® Agar - Dehydrated base medium                  | 5 kg      | 5 years    | 455                                 | 2-30°C     | -                                                                                          | -            | -                    |
| <b>BK216HA</b>  | EASY Staph® Agar - Dehydrated base medium                  | 500 g     | 5 years    | 455                                 | 2-30°C     | -                                                                                          | -            | -                    |
| <b>BK219HA</b>  | Alcaline Saline Peptone water                              | 500 g     | 5 years    | 455                                 | 2-30°C     | 6 months at 2-25°C                                                                         | 15 min 121°C | -                    |
| <b>BK221HA</b>  | TSI Agar (ISO 6579-1)                                      | 500 g     | 5 years    | 455                                 | 2-30°C     | in tubes not inclined (*) : 180 days at 2-8 °C.<br>in inclined tubes (*) : 8 days at 2-8 ° | -            | -                    |
| <b>BK223HA</b>  | Bismuth-Sulfite (ISO 6579-1) Agar                          | 500 g     | 5 years    | 455                                 | 2-30°C     | 4 days at 2-8°C, shielded from light                                                       | 0            | -                    |
| <b>BK224HA</b>  | TSYE Agar (ISO 11290-1)                                    | 500 g     | 5 years    | 455                                 | 2-30°C     | 180 days at 2-25°C                                                                         | -            | -                    |
| <b>BK225GC</b>  | <i>Salmonella</i> Enrichment double strength buffered      | 5 kg      | 5 years    | 455                                 | 2-30°C     | 2-25°C                                                                                     | 0            | -                    |
| <b>BK225HA</b>  | Salmonella Enrichment double strength buffered             | 500 g     | 5 years    | 455                                 | 2-30°C     | 2-25°C                                                                                     | 0            | -                    |
| <b>BK227GC</b>  | SYMPHONY® Agar                                             | 5 kg      | 5 years    | 455                                 | 2-30°C     | Fl. Tb. 180 days at 2-25°C; BP 30 days at 2-8 °C.                                          | 0            | -                    |
| <b>BK227HA</b>  | SYMPHONY® Agar                                             | 500 g     | 5 years    | 455                                 | 2-30°C     | Fl. Tb. 180 days at 2-25°C; BP 30 days at 2-8 °C.                                          | 0            | -                    |



| Designation ENG | Packaging                                                               | Packaging                    | Shelf life | Minimal Remaining shelf life (days) | Storing T°                 | Storage after preparation             | Autoclave    | Supplements |
|-----------------|-------------------------------------------------------------------------|------------------------------|------------|-------------------------------------|----------------------------|---------------------------------------|--------------|-------------|
| <b>BK232GC</b>  | Chromoclin Agar Urine                                                   | 5 kg                         |            |                                     | 2-30°C                     |                                       |              |             |
| <b>BK232HA</b>  | Chromoclin Agar Urine                                                   | 500 g                        |            |                                     | 2-30°C                     |                                       |              |             |
| <b>BK239HA</b>  | Iron Sulfite Agar (ISA) (ISO 15213-1)                                   | 500 g                        | 5 years    | 455                                 | 2-30°C                     | 180 days at 2-25°C                    | 15 min 121°C | -           |
| <b>BK241HA</b>  | SIM Agar                                                                | 500 g                        | 3 years    | 455                                 | 2-20°C                     |                                       |              |             |
| <b>BK242GC</b>  | Buffered-peptone-water-HALAL (BPW HALAL) ISO                            | 5 kg                         | 5 years    | 455                                 | 2-30°C                     |                                       |              |             |
| <b>BK242HA</b>  | Buffered-peptone-water-HALAL (BPW HALAL) ISO                            | 500 g                        | 5 years    | 455                                 | 2-30°C                     |                                       |              |             |
| <b>BK243GC</b>  | Trypto casein soy broth HALAL(TSB HLAL) ISO/EP/USP/JP                   | 5 kg                         | 5 years    | 455                                 | 2-30°C                     |                                       |              |             |
| <b>BK243HA</b>  | Trypto casein soy broth HALAL(TSB HLAL) ISO/EP/USP/JP                   | 5 kg                         | 5 years    | 455                                 | 2-30°C                     |                                       |              |             |
| <b>BK244GC</b>  | Trypto-caseine soy agar HALAL (TSA HALAL) ISO/EP/USP/JP                 | 500 g                        | 5 years    | 455                                 | 2-30°C                     |                                       |              |             |
| <b>BK244HA</b>  | Trypto-caseine soy agar HALAL (TSA HALAL) ISO/EP/USP/JP                 | 5 kg                         | 5 years    | 455                                 | 2-30°C                     |                                       |              |             |
| <b>BK245GC</b>  | EP/USP/JP                                                               | 500 g                        | 5 years    | 455                                 | 2-30°C                     |                                       |              |             |
| <b>BK245HA</b>  | EP/USP/JP                                                               | 5 kg                         | 5 years    | 455                                 | 2-30°C                     |                                       |              |             |
| <b>BK248HA</b>  | M17 S Agar                                                              | 500 g                        | 5 years    | 455                                 | 2-20°C                     |                                       |              |             |
| <b>BK250HA</b>  | Oxford Agar - Dehydrated base medium                                    | 500 g                        | 5 years    | 455                                 | 2-30°C                     | Base 6 months, plate 30 days at 2-8°C |              |             |
| <b>BM00608</b>  | Eugon LT 100 Broth - Ready-to-use medium                                | 50 × 9 mL in tubes           | 365 days   | 90                                  | 2-25°C                     | -                                     | -            | -           |
| <b>BM00808</b>  | Tryptone-salt Broth (Maximum Recovery Diluent) - Ready-to-use medium    | 50 × 9 mL in tubes           | 548 days   | 90                                  | 2-25°C                     | -                                     | -            | -           |
| <b>BM00908</b>  | Trypto-casein Soy Broth (TSB) - Ready-to-use medium                     | 10 × 100 mL in vials         | 548 days   | 90                                  | 2-25°C                     | -                                     | -            | -           |
| <b>BM01008</b>  | Buffered peptone water (25.5 g/L) - Ready-to-use medium                 | 10 × 225 mL in vials         | 365 days   | 90                                  | 2-25°C                     | -                                     | -            | -           |
| <b>BM01108</b>  | Brilliant Green Bile Broth (BGBB) single strength - Ready-to-use medium | 50 × 10 mL with Durham tubes | 365 days   | 90                                  | 2-25°C                     | -                                     | -            | -           |
| <b>BM01308</b>  | FRASER Broth - Ready-to-use medium                                      | 50 × 10 mL in tubes          | 365 days   | 90                                  | 2-8°C, shielded from light | -                                     | -            | -           |
| <b>BM01508</b>  | Plate count Agar (PCA) - Ready-to-melt medium                           | 10 × 100 mL in vials         | 548 days   | 90                                  | 2-25°C                     | -                                     | -            | -           |
| <b>BM01608</b>  | Half-FRASER Broth - Ready-to-use medium                                 | 10 × 225 mL in vials         | 365 days   | 90                                  | 2-8°C, shielded from light | -                                     | -            | -           |

| Designation ENG | Packaging                                                         | Packaging            | Shelf life | Minimal Remaining shelf life (days) | Storing T° | Storage after preparation | Autoclave | Supplements |
|-----------------|-------------------------------------------------------------------|----------------------|------------|-------------------------------------|------------|---------------------------|-----------|-------------|
| <b>BM01708</b>  | Trypto-casein Soy Agar (TSA) - Ready-to-melt medium               | 10 × 100 mL in vials | 548 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM01808</b>  | BAIRD PARKER Agar with Egg Yolk and Tellurite - Pre-poured medium | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C      | -                         | -         | -           |
| <b>BM02008</b>  | PALCAM Agar - Pre-poured medium                                   | 20 plates Ø 90 mm    | 150 days   | 60                                  | 2-8°C      | -                         | -         | -           |
| <b>BM02108</b>  | Chloramphenicol glucose Agar - Ready-to-melt medium               | 10 × 100 mL in vials | 365 days   | 90                                  | 2-8°C      | -                         | -         | -           |
| <b>BM02208</b>  | Oxytetracycline glucose Agar (OGA) - Ready-to-melt base medium    | 10 x 110 mL in vials | 365 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM03008</b>  | Trypto-casein Soy Broth (TSB) - Ready-to-use medium               | 50 × 10 mL in tubes  | 548 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM03308</b>  | Plate count Agar (PCA) - Ready-to-melt medium                     | 10 × 200 mL in vials | 548 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM03408</b>  | VRBL Agar - Ready-to-melt medium                                  | 10 × 100 mL in vials | 365 days   | 90                                  | 2-8°C      | -                         | -         | -           |
| <b>BM03508</b>  | VRBL Agar - Ready-to-melt medium                                  | 10 × 200 mL in vials | 365 days   | 90                                  | 2-8°C      | -                         | -         | -           |
| <b>BM03608</b>  | XLT4 Agar - Pre-poured medium                                     | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C      | -                         | -         | -           |
| <b>BM03808</b>  | Bacillus cereus (acc. to Mossel) - Pre-poured medium              | 20 plates Ø 90 mm    | 150 days   | 60                                  | 2-8°C      | -                         | -         | -           |
| <b>BM03908</b>  | TSC Agar - Ready-to-use base medium                               | 50 × 20 mL in tubes  | 365 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM04308</b>  | Eugon LT 100 Broth - Ready-to-use medium                          | 10 × 100 mL in vials | 365 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM04908</b>  | Trypto-casein Soy Agar (TSA) - Ready-to-melt medium               | 10 × 200 mL in vials | 365 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM05008</b>  | Trypto-casein Soy Agar (TSA) - Pre-poured medium                  | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C      | -                         | -         | -           |
| <b>BM05308</b>  | Sabouraud Dextrose Agar (SDA) - Ready-to-melt medium              | 10 x 200 mL in vilas | 365 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM05608</b>  | Buffered peptone water (25.5 g/L) - Ready-to-use medium           | 50 × 9 mL in tubes   | 365 days   | 90                                  | 2-25°C     | -                         | -         | -           |
| <b>BM05708</b>  | Buffered peptone water (25.5 g/L) - Ready-to-use medium           | 10 × 90 mL in vials  | 365 days   | 90                                  | 2-25°C     | -                         | -         | -           |

| Designation ENG | Packaging                                                      | Packaging            | Shelf life | Minimal Remaining shelf life (days) | Storing T°                 | Storage after preparation | Autoclave | Supplements |
|-----------------|----------------------------------------------------------------|----------------------|------------|-------------------------------------|----------------------------|---------------------------|-----------|-------------|
| <b>BM06608</b>  | COMPASS® <i>Salmonella</i> Agar - Pre-poured medium            | 20 plates Ø 90 mm    | 90 days    | 30                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM06708</b>  | BAIRD-PARKER RPF Agar - Pre-poured medium                      | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM06808</b>  | Yeast extract Agar - Ready-to-melt medium                      | 10 × 200 mL in vials | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM06908</b>  | TBX Agar Ready-to-melt medium                                  | 10 × 100 mL in vials | 270 days   | 90                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM07108</b>  | GVPC Agar for <i>Legionella</i> - Pre-poured medium            | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM07208</b>  | BCYE Agar for <i>Legionella</i> - Pre-poured medium            | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM07308</b>  | BCYE Agar for <i>Legionella</i> - Pre-poured medium            | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM07408</b>  | RAPPAPORT-VASSILIADIS Soja (RVS) Broth - Ready-to-use medium   | 50 × 10 mL in tubes  | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM07508</b>  | VRBG Agar - Ready-to-melt medium                               | 10 × 200 mL in vials | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM07608</b>  | Tryptophan Broth - Ready-to-use medium                         | 50 × 3 mL in tubes   | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM07708</b>  | TSC Agar - Ready-to-melt base medium                           | 10 × 200 mL in vials | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM07808</b>  | MKTTn Broth - Ready-to-use medium                              | 50 × 10 mL in tubes  | 548 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM07908</b>  | Chloramphenicol glucose Agar - Ready-to-melt medium            | 10 × 200 mL in vials | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM08208</b>  | Thioglycollate medium with resazurin - Ready-to-use medium     | 50 × 10 mL in tubes  | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM08608</b>  | Plate count Agar with skimmed milk - Ready-to-melt medium      | 10 × 200 mL in vials | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM08708</b>  | XLD Agar (ISO 6579-1) - Pre-poured medium                      | 20 plates Ø 90 mm    | 90 days    | 30                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM08808</b>  | Synthetic sea salt with Bromophenol blue - Ready-to-use medium | 50 × 18 mL in tubes  | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |

| Designation ENG | Packaging                                                                                    | Packaging                    | Shelf life | Minimal Remaining shelf life (days) | Storing T°                 | Storage after preparation | Autoclave | Supplements |
|-----------------|----------------------------------------------------------------------------------------------|------------------------------|------------|-------------------------------------|----------------------------|---------------------------|-----------|-------------|
| <b>BM08908</b>  | MRS Agar pH 5.7 - Ready-to-melt medium                                                       | 10 × 200 mL in vials         | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM09108</b>  | BAIRD PARKER Agar with Egg Yolk and Tellurite - Pre-poured medium                            | 120 plates Ø 90 mm           | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM09308</b>  | Tergitol 7 Agar - Pre-poured medium                                                          | 120 plates Ø 55 mm           | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM09408</b>  | SLANETZ and BARTLEY Agar - Pre-poured medium                                                 | 120 plates Ø 55 mm           | 150 days   | 60                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM09508</b>  | Acetamide Broth - Ready-to-use                                                               | 7 × 5 mL in tubes            | 365 days   | 90                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM09608</b>  | CFC Agar - Ready-to-melt base medium                                                         | 10 × 200 mL in vials         | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM09708</b>  | Laurylsulfate-Tryptose single strength Broth - Ready-to-use medium                           | 50 × 10 mL with Durham tubes | 365 days   | 90                                  | 15-25°C                    | -                         | -         | -           |
| <b>BM09808</b>  | Laurylsulfate-Tryptose double strength Broth - Ready-to-use medium                           | 50 × 10 mL with Durham tubes | 365 days   | 90                                  | 15-25°C                    | -                         | -         | -           |
| <b>BM09908</b>  | Glucose Agar - Ready-to-use medium                                                           | 50 × 10 mL in tubes          | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM10408</b>  | BEA Agar - Ready-to-melt medium                                                              | 10 × 100 mL in vials         | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM10508</b>  | KING B Agar - Ready-to-melt medium                                                           | 7 × 7 mL in tubes            | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM10808</b>  | TSYE Agar - Ready-to-melt medium                                                             | 50 × 18 mL in tubes          | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM10908</b>  | Dichloran-Glycerol (DG18) Agar - Ready-to-melt medium                                        | 10 × 100 mL in vials         | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM11008</b>  | GIOLITTI and CANTONI with Tween®80 single strength concentration Broth - Ready-to-use medium | 50 × 10 mL in tubes          | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM11108</b>  | GIOLITTI and CANTONI with Tween®80 double strength concentration Broth - Ready-to-use medium | 50 × 10 mL in tubes          | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |

| Designation ENG | Packaging                                                            | Packaging            | Shelf life | Minimal Remaining shelf life (days) | Storing T°                 | Storage after preparation | Autoclave | Supplements |
|-----------------|----------------------------------------------------------------------|----------------------|------------|-------------------------------------|----------------------------|---------------------------|-----------|-------------|
| <b>BM11408</b>  | Tryptone-salt Broth (Maximum Recovery Diluent) - Ready-to-use medium | 10 × 90 mL in vials  | 548 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM11508</b>  | Sterile distilled water                                              | 50 × 18 mL in tubes  | 548 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM11608</b>  | COMPASS® Enterococcus Agar - Ready-to-melt medium                    | 10 × 100 mL in vials | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM11808</b>  | Nutrient Agar 2% - Ready-to-melt medium                              | 50 × 18 mL in tubes  | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM12208</b>  | Sugar-Free Agar (SFA) - Ready-to-melt medium                         | 10 × 100 mL in vials | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM12308</b>  | COMPASS® <i>Listeria</i> Agar - Pre-poured medium                    | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM12408</b>  | COMPASS® <i>Listeria</i> Agar - Pre-poured medium                    | 120 plates Ø 90 mm   | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM12508</b>  | Nutrient Agar 2.5% - Ready-to-melt medium                            | 50 × 18 mL in tubes  | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM12608</b>  | COMPASS® <i>Bacillus cereus</i> Agar - Pre-poured medium             | 20 plates Ø 90 mm    | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM12708</b>  | MSRV Agar (ISO 6579-1) - Ready-to-melt medium                        | 200 mL in vials      | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM13008</b>  | COMPASS® <i>Bacillus cereus</i> Agar (Base) - Ready-to-melt medium   | 100 mL in vials      | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM13108</b>  | Buffered peptone water (25.5 g/L) - Ready-to-use medium              | L in flexible        | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM13208</b>  | Buffered peptone water (25.5 g/L) - Ready-to-use medium              | L in flexible        | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM13308</b>  | Half-FRASER Broth - Ready-to-use medium                              | L in flexible        | 365 days   | 90                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM13408</b>  | Half-FRASER Broth - Ready-to-use medium                              | L in flexible        | 365 days   | 90                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM13508</b>  | Tryptone-salt Broth (Maximum Recovery Diluent) - Ready-to-use medium | L in flexible        | 548 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM13608</b>  | <i>Salmonella</i> Enrichment - Ready-to-use medium                   | 225 mL in vials      | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM13708</b>  | <i>Salmonella</i> Enrichment - Ready-to-use medium                   | L in flexible        | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM13808</b>  | SESAME <i>Salmonella</i> Detection - Ready-to-melt medium            | 200 mL in vials      | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM13908</b>  | CONFIRM' L. mono Agar - Pre-poured medium                            | plates Ø 90 mm       | 120 days   | 30                                  | 2-8°C                      | -                         | -         | -           |

| Designation ENG | Packaging                                                                | Packaging     | Shelf life | Minimal Remaining shelf life (days) | Storing T°                 | Storage after preparation | Autoclave | Supplements |
|-----------------|--------------------------------------------------------------------------|---------------|------------|-------------------------------------|----------------------------|---------------------------|-----------|-------------|
| <b>BM14208</b>  | Dichloran Rose Bengal Chloramphenicol (DRBC) Agar - Ready-to-melt medium | x 200 mL in v | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM14408</b>  | <i>Salmonella</i> Enrichment - Ready-to-use medium                       | L in flexible | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM14508</b>  | CN Agar for Pseudomonas - Pre-poured medium                              | plates Ø 55   | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM14608</b>  | SLANETZ and BARTLEY Agar - Pre-poured medium                             | plates Ø 55   | 150 days   | 60                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM14708</b>  | Tergitol 7 Agar - Pre-poured medium                                      | plates Ø 55   | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM15008</b>  | SESAME <i>Salmonella</i> Detection - Pre-poured medium                   | plates Ø 90   | 150 days   | 30                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM15308</b>  | COMPASS® cc Agar - Pre-poured medium                                     | plates Ø 55   | 120 days   | 40                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM15408</b>  | <i>Cronobacter</i> Chromogenic Isolation (CCI) agar - Pre-poured medium  | plates Ø 90   | 90 days    | 30                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM15508</b>  | <i>Cronobacter</i> Screening Broth (CSB) - Ready-to-use medium           | x 10 mL in tu | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM15608</b>  | Penicillin Pimaricine Agar (PPA) - Pre-poured medium                     | plates Ø 90   | 90 days    | 30                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM15708</b>  | COMPASS® Enterococcus Agar - Pre-poured medium                           | plates Ø 55   | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM15908</b>  | BAIRD-PARKER RPF Agar - Pre-poured medium                                | plates Ø 55   | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM16008</b>  | IRIS <i>Salmonella</i> ® Agar - Pre-poured medium                        | plates Ø 90   | 90 days    | 30                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM16108</b>  | IRIS <i>Salmonella</i> ® Agar - Pre-poured medium                        | plates Ø 90   | 90 days    | 30                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM16208</b>  | CONFIRM' L. mono Broth                                                   | x 1 mL in vi  | 270 days   | 90                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM16308</b>  | <i>Salmonella</i> Enrichment with Tween®80 - Ready-to-use                | L in flexible | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM16708</b>  | RHAPSODY Agar® - Pre-poured medium                                       | plates Ø 90   | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM16808</b>  | BCP Dextrose Agar - Ready-to-melt medium                                 | x 200 mL in v | 365 days   | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM16908</b>  | Meat liver glucose agar with 2g/L yeast extract                          | x 200 mL in v | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM17108</b>  | TBX Agar Ready-to-melt medium                                            | x 200 mL in v | 270 days   | 90                                  | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM17208</b>  | Sabouraud Dextrose Chloramphenicol Agar (SDCA) - Ready-to-melt medium    | x 200 mL in v | 365 days   | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM17308</b>  | Sabouraud Dextrose Agar (SDA) - Pre-poured medium                        | plates Ø 90   | 120 days   | 40                                  | 2-8°C                      | -                         | -         | -           |



| Designation ENG | Packaging                                                             | Packaging                 | Shelf life | Minimal<br>Remaining<br>shelf life<br>(days) | Storing T°                   | Storage after preparation | Autoclave | Supplements |
|-----------------|-----------------------------------------------------------------------|---------------------------|------------|----------------------------------------------|------------------------------|---------------------------|-----------|-------------|
| <b>BM17608</b>  | Sabouraud Dextrose Chloramphenicol Agar (SDCA) - Ready-to-melt medium | 20 plates Ø 90 mm         | 120 days   | 40                                           | 2-8°C                        | -                         | -         | -           |
| <b>BM17808</b>  | Trypto-casein Soy Agar (TSA) - Pre-poured medium                      | 120 plates Ø 90 mm        | 120 days   | 40                                           | 2-8°C                        | -                         | -         | -           |
| <b>BM17908</b>  | Trypto-casein Soy Broth (TSB) - Ready-to-use medium                   | 10 x 90 mL in vials       | 548 days   | 90                                           | 2-25 °C, shielded from light | -                         | -         | -           |
| <b>BM18008</b>  | MacCONKEY Agar - Ready-to-melt medium                                 | 10 x 100 mL in vials      | 365 days   | 90                                           | 2-8°C                        | -                         | -         | -           |
| <b>BM18208</b>  | Chromogenic Coliform Agar (CCA - ISO 9308-1)                          | 20 plates Ø 55 mm         | 120 days   | 40                                           | 2-8°C, shielded from light   | -                         | -         | -           |
| <b>BM18308</b>  | R2A - Ready-to-melt medium                                            | 10 x 100 mL in vials      | 365 days   | 90                                           | 2-25°C                       | -                         | -         | -           |
| <b>BM18408</b>  | Cetrimide Agar - Ready-to-melt medium                                 | 10 x 100 mL in vials      | 270 days   | 90                                           | 2-8°C                        | -                         | -         | -           |
| <b>BM18608</b>  | Potato Dextrose Agar (PDA) - Ready-to-melt medium                     | 10 x 200 mL in vials      | 365 days   | 90                                           | 2-25°C                       | -                         | -         | -           |
| <b>BM18708</b>  | EASY <i>Staph</i> ® Agar - pre-poured                                 | 20 plates Ø 90 mm         | 120 days   | 40                                           | 2-8°C                        | -                         | -         | -           |
| <b>BM18808</b>  | Half-FRASER Broth - Ready-to-use medium                               | 40 x 5 L in flexible bags | 365 days   | 90                                           | 2-8°C, shielded from light   | -                         | -         | -           |
| <b>BM19008</b>  | EASY <i>Staph</i> ® Agar - pre-poured                                 | 120 plates Ø 90 mm        | 120 days   | 40                                           | 2-8°C                        | -                         | -         | -           |
| <b>BM19108</b>  | SYMPHONY® Agar - Ready-to-melt medium                                 | 10 x 200 mL in vials      | 365 days   | 90                                           | 2-8°C                        | -                         | -         | -           |
| <b>BM19208</b>  | Lactose-sulfite (LS) Broth - Ready-to-use medium                      | 7 x 9 mL in tubes         | 120 days   | 40                                           | 2-8°C, shielded from light   | -                         | -         | -           |
| <b>BM19308</b>  | Dipotassium hydrogenphosphate solution - Ready-to-use medium          | 2 x 5 L in flexible bags  | 18 months  | 90                                           | 2-25°C                       | -                         | -         | -           |
| <b>BM19408</b>  | Sterile distilled water                                               | 10 x 90 mL in vials       | 548 jours  | 90                                           | 2-25°C                       | -                         | -         | -           |
| <b>BM19608</b>  | CN Agar for Pseudomonas - Pre-poured medium                           | 120 plates Ø 55 mm        | 120 days   | 40                                           | 2-8°C                        | -                         | -         | -           |

| Designation ENG | Packaging                                                                | Packaging                 | Shelf life | Minimal Remaining shelf life (days) | Storing T°                   | Storage after preparation | Autoclave | Supplements |
|-----------------|--------------------------------------------------------------------------|---------------------------|------------|-------------------------------------|------------------------------|---------------------------|-----------|-------------|
| <b>BM19708</b>  | OF glucose Medium                                                        | 50 × 10 mL in tubes       | 365 days   | 90                                  | 2-8°C                        | -                         | -         | -           |
| <b>BM19808</b>  | <i>Salmonella</i> Enrichment with Tween®80 - Ready-to-use                | 2 × 5 L in flexible bags  | 365 days   | 90                                  | 2-25°C                       | -                         | -         | -           |
| <b>BM19908</b>  | Bacillus cereus Agar (acc. to Mossel) - Ready-to-use                     | 120 plates Ø 90 mm        | 150 days   | 60                                  | 2-8°C                        | -                         | -         | -           |
| <b>BM20008</b>  | <i>Salmonella</i> Enrichment double strength buffered - Ready-to-use     | 2 × 5 L in flexible bags  | 365 days   | 90                                  | 2-25°C                       | -                         | -         | -           |
| <b>BM20108</b>  | <i>Salmonella</i> Enrichment double strength buffered - Ready-to-use     | 10 x 225 mL in vials      | 365 days   | 90                                  | 2-25°C                       | -                         | -         | -           |
| <b>BM20208</b>  | SYMPHONY® Agar - Pre-poured medium                                       | 20 plates Ø 90 mm         | 120 days   | 40                                  | 2-8°C                        | -                         | -         | -           |
| <b>BM20608</b>  | PCA + TTC + Neutralizers - Contact Plate                                 | 20 plates Ø 65 mm         | 120 jours  | 30                                  | 2-8°C                        | -                         | -         | -           |
| <b>BM20609</b>  | CONTACT PCA + TTC + NEUTRALIZERS                                         | 30 plates                 | 6 months   |                                     | 2-25°C                       |                           |           |             |
| <b>BM20708</b>  | VRBG + Neutralizers - Contact Plate                                      | 20 plates Ø 65 mm         | 120 jours  | 30                                  | 2-8°C                        | -                         | -         | -           |
| <b>BM20709</b>  | CONTACT VRBG + NEUTRALIZERS                                              | 30 plates                 | 7 monthes  |                                     | 2-25°C                       |                           |           |             |
| <b>BM20908</b>  | VRBL + Neutralizers - Contact Plate                                      | 20 plates Ø 65 mm         | 90 days    | 30                                  | 2-8°C                        | -                         | -         | -           |
| <b>BM20909</b>  | CONTACT VRBL + NEUTRALIZERS                                              | 30 plates                 | 7 monthes  |                                     | 2-25°C                       |                           |           |             |
| <b>BM21008</b>  | Sabouraud Dextrose Chloramphenicol (SDCA) + Neutralizers - Contact plate | 20 plates Ø 65 mm         | 120 days   | 30                                  | 2-8°C (unopened box at 25°C) | -                         | -         | -           |
| <b>BM21108</b>  | BAIRD-PARKER + Neutralizers - Contact plate                              | 20 plates Ø 65 mm         | 120days    | 30                                  | 2-8°C (unopened box at 25°C) | -                         | -         | -           |
| <b>BM21109</b>  | CONTACT BAIRD-PARKER + NEUTRALIZERS                                      | 30 plates                 | 7 monthes  |                                     | 2-14°C                       |                           |           |             |
| <b>BM21208</b>  | Half Fraser + Tween® 80 (10 g/L)                                         | 1 L in flexible bag       | 1 year     | 90                                  | 2-8°C, shielded from light   | -                         | -         | -           |
| <b>BM21308</b>  | COMPASS® Bacillus Plus Agar                                              | 20 plates Ø 90 mm         | 120days    | 40                                  | 2-8°C                        | -                         | -         | -           |
| <b>BM21408</b>  | Buffered peptone water (25.5 g/L) - Ready-to-use medium                  | 40 × 5 L in flexible bags | 1 year     | 90                                  | 2-25°C                       | -                         | -         | -           |

| Designation ENG | Packaging                                                          | Packaging                           | Shelf life              | Minimal Remaining shelf life (days) | Storing T°                 | Storage after preparation | Autoclave | Supplements |
|-----------------|--------------------------------------------------------------------|-------------------------------------|-------------------------|-------------------------------------|----------------------------|---------------------------|-----------|-------------|
| <b>BM21708</b>  | Plate count Agar (PCA) - Pre-poured medium                         | 20 plates Ø 90 mm                   | 120 days                | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM21808</b>  | Plate count Agar (PCA) - Pre-poured medium                         | 120 plates Ø 90 mm                  | 120 days                | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM21908</b>  | Eugon LT Sup- broth                                                | 50 x 9 mL in tubes                  | 365 days                | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM22208</b>  | Eugon LT Sup- broth                                                | 10 x 100 mL in vials                | 365 days                | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM22408</b>  | Eugon LT Sup- broth                                                | 50 x 9 mL in tubes with glass beads | 365 days                | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM22708</b>  | GVPC Agar for <i>Legionella</i> - Pre-poured medium                | 120 plates Ø 90 mm                  | 120 days                | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM22808</b>  | <i>Salmonella</i> Enrichment with Tween®80 - Ready-to-use          | 10 x 225 mL in vials                | 1 year                  | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM22908</b>  | Eugon LT Sup- broth                                                | 5 x 5 L in flexible bags            | 365 days                | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM23008</b>  | COMPASS® <i>Salmonella</i> Agar pre-poured medium                  | 120 plates Ø 90 mm                  | 90 days                 | 30                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM23208</b>  | PERFECTO Duo Agar - ready to melt base medium                      | x200ml in vial                      | Not created yet         | pas encore créée                    | 2-8°C, shielded from light | -                         | -         | -           |
| <b>BM23308</b>  | <i>Salmonella</i> enrichment broth (BPW conc.) 9X                  | 5L in flexible                      | stability test in going | 90                                  | 15-25 °C                   | -                         | -         | -           |
| <b>BM23508</b>  | Iron Sulfite Agar (ISA) (ISO 15213-1) Ready-to-melt medium         | 10 x 200 mL in vials                | 1 year                  | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM23608</b>  | COMPASS® <i>Bacillus cereus</i> Agar (base) - Ready-to-melt medium | 10 x 200 mL in vials                | 365 days                | 90                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM23708</b>  | <i>Salmonella</i> Enrichment - Ready-to-use medium                 | 40 x 5 L in flexible bags           | 365 days                | 90                                  | 2-25°C                     | -                         | -         | -           |
| <b>BM23908</b>  | MRS pre-poured medium                                              | 120 plates Ø 90 mm                  | -                       | -                                   | 2-8°C                      | -                         | -         | -           |
| <b>BM24108</b>  | SIM Agar ready-to-use                                              | x 10 mL in tube                     | 6 months                | 60                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM24308</b>  | SYMPHONY® Agar - Pre-poured medium                                 | 120 plates Ø 90 mm                  | 120 days                | 40                                  | 2-8°C                      | -                         | -         | -           |
| <b>BM24409</b>  | CONTACT SDCA + NEUTRALIZERS                                        | 30 plates                           | 7 months                |                                     | 2-25°C                     |                           |           |             |

| Designation ENG | Packaging                                  | Packaging                         | Shelf life          | Minimal<br>Remaining<br>shelf life<br>(days) | Storing T°                    | Storage after preparation             | Autoclave | Supplements |
|-----------------|--------------------------------------------|-----------------------------------|---------------------|----------------------------------------------|-------------------------------|---------------------------------------|-----------|-------------|
| <b>BR00108</b>  | Ringer's solution (1/4 strength)           | 100 tablets                       | 3 years / 6 months* | 365                                          | 2-30°C                        | 180 days at 2-25°C                    | -         | -           |
| <b>BR00208</b>  | Lyophilized coagulase rabbit plasma        | 10 vials<br>pack for 20 reactions | 2 years / 30 days*  | 240                                          | 2-8°C                         | 30 days at 2-8°C                      | 0         | -           |
| <b>BR00308</b>  | Synthetic sea salt                         | 100 g                             | 5 years             | 365                                          | 2-30°C                        | 180 days at 2-25°C                    | -         | -           |
| <b>BS00408</b>  | PALCAM selective supplement                | 10 vials<br>pack qsf 500 ml       | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 days at 2-8°C, shielded from light | 0         | -           |
| <b>BS00608</b>  | D-cycloserine 200 mg selective supplement  | 10 vials<br>pack qsf 500 ml       | 3 years / 20 days*  | 365                                          | 2-8°C                         | 20 days at 2-8°C                      | 0         | -           |
| <b>BS00708</b>  | Polymyxin B selective supplement           | 10 vials<br>pack                  | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 days at 2-8°C                      | 0         | -           |
| <b>BS00808</b>  | Oxytetracycline 50 mg selective supplement | 10 vials<br>pack qsf 550 ml       | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 days at 2-8°C                      | 0         | -           |
| <b>BS00908</b>  | Gentamicin 25 mg selective supplement      | 10 vials<br>pack                  | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 days at 2-8°C                      | 0         | -           |
| <b>BS02108</b>  | Chloramphenicol 50 mg selective supplement | 10 vials<br>pack                  | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 at 2-8 °C                          | 0         | -           |
| <b>BS02208</b>  | CFC selective supplement                   | 10 vials<br>pack qsf 500 ml       | 2 years / 30 days*  | 365                                          | 2-8°C                         | 30 at 2-8 °C                          | 0         | -           |
| <b>BS02408</b>  | MUG 50 mg supplement                       | 10 vials<br>pack                  | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 at 2-8 °C                          | 0         | -           |
| <b>BS02608</b>  | TTC supplement 12.5 mg                     | 10 vials<br>pack                  | 3 years / 30 days*  | 365                                          | 2-8°C,<br>shielded from light | 30 at 2-8°C shielded from light       | 0         | -           |
| <b>BS02708</b>  | TTC 50 mg supplement                       | 10 vials                          | 3 years / 30 days*  | 365                                          | 2-8°C,<br>shielded from light | 30 at 2-8°C shielded from light       | 0         | -           |
| <b>BS02808</b>  | SULFAMETHAZINE 25 mg SELECTIVE SUPPLEMENT  | ls pack qsf 5                     | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 at 2-8 °C                          | 0         | -           |
| <b>BS03008</b>  | HALF-FRASER SELECTIVE SUPPLEMENT           | ls pack qsf 5                     | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 at 2-8°C shielded from light       | 0         | -           |
| <b>BS03108</b>  | FRASER SELECTIVE SUPPLEMENT                | ls pack qsf 5                     | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 at 2-8°C shielded from light       | 0         | -           |
| <b>BS03208</b>  | HALF-FRASER SELECTIVE SUPPLEMENT           | ls pack qsf 2                     | 3 years / 30 days*  | 365                                          | 2-8°C                         | 30 at 2-8°C shielded from light       | 0         | -           |

| Designation ENG | Packaging                                                           | Packaging                        | Shelf life         | Minimal Remaining shelf life (days) | Storing T°                 | Storage after preparation       | Autoclave | Supplements |
|-----------------|---------------------------------------------------------------------|----------------------------------|--------------------|-------------------------------------|----------------------------|---------------------------------|-----------|-------------|
| <b>BS03308</b>  | Novobiocin 10 mg selective supplement                               | 10 vials                         | 3 years / 30 days* | 365                                 | 2-8°C                      | 30 at 2-8 °C                    | 0         | -           |
| <b>BS03408</b>  | Rabbit Plasma Fibrinogen (RPF) supplement                           | 8 vials pack (1 vial qsf 100 mL) | 3 years / 8 days*  | 365                                 | 2-8°C                      | 8 at 2-8°C                      | 0         | -           |
| <b>BS03608</b>  | Sterile Egg yolk emulsion with Potassium tellurite                  | 1 vial qsf 900 mL                | 2 years / 30 days* | 240                                 | 2-8°C, shielded from light | 30 at 2-8 °C                    | 0         | -           |
| <b>BS03708</b>  | Cefixime-tellurite supplement for CT-SMAC Agar                      | 10 vials pack qsf 500 mL         | 2 years / 30 days* | 365                                 | 2-8°C                      | 30 at 2-8°C                     | 0         | -           |
| <b>BS03808</b>  | Rabbit Plasma Fibrinogen (RPF) supplement                           | 1 vial qsf 500 mL                | 3 years / 8 days*  | 365                                 | 2-8°C                      | 8 at 2-8°C                      | 0         | -           |
| <b>BS04908</b>  | PALCAM selective supplement                                         | 8 vials pack qsf 2.5 L           | 3 years / 30 days* | 365                                 | 2-8°C                      | 30 at 2-8°C shielded from light | 0         | -           |
| <b>BS05508</b>  | Egg Yolk Emulsion with Polymyxin B                                  | 10 vials pack qsf 500 mL         | 2 years            | 240                                 | 2-8°C                      | -                               | -         | -           |
| <b>BS05608</b>  | Novobiocin 40 mg selective supplement                               | 10 vials pack                    | 3 years / 30 days* | 365                                 | 2-8°C                      | 30 at 2-8 °C                    | 0         | -           |
| <b>BS05908</b>  | Sterile ferric ammonium citrate 5% supplement                       | 10 × 90 mL in vials              | 2 years            | 240                                 | 2-25°C                     | -                               | -         | -           |
| <b>BS06008</b>  | Sterile Egg Yolk Tellurite Enrichment                               | 10 × 50 mL in vials              | 2 years            | 240                                 | 2-8°C, shielded from light | -                               | -         | -           |
| <b>BS06208</b>  | Sterile ferric ammonium citrate 5% supplement                       | 7 × 10 mL in tubes               | 2 years / 30 days* | 240                                 | 2-25°C                     | -                               | -         | -           |
| <b>BS06608</b>  | Egg Yolk Enrichment                                                 | 10 × 50 mL in vials              | 2 years            | 240                                 | 2-8°C                      | -                               | -         | -           |
| <b>BS06908</b>  | COMPASS® Bacillus cereus - Selective supplement                     | 10 vials pack qsf 100 mL         | 2 years / 30 days* | 240                                 | 2-8°C                      | 30 °C                           | 0         | -           |
| <b>BS07008</b>  | COMPASS® <i>Listeria</i> Agar - Enrichment supplement (liquid)      | 8 vials pack qsf 1 L             | 1 year             | 90                                  | 2-25°C                     | -                               | -         | -           |
| <b>BS07108</b>  | COMPASS® <i>Listeria</i> Agar - Selective supplement (lyophilisate) | 8 vials pack qsf 1 L             | 2 years / 15 days* | 240                                 | 2-8°C                      | 15 at 2-8°C shielded from light | 0         | -           |
| <b>BS07708</b>  | IRIS <i>Salmonella</i> ® supplement qsf 225 mL - Tablets            | 120 tablets                      | 42 months          | 180                                 | 2-8°C                      | -                               | -         | -           |
| <b>BS07808</b>  | IRIS <i>Salmonella</i> ® supplement - Liquid                        | 10 × 50 mL in vials              | 9 months           | 90                                  | 2-8°C                      | -                               | -         | -           |

| Designation ENG | Packaging                                               | Packaging                    | Shelf life         | Minimal Remaining shelf life (days) | Storing T°                 | Storage after preparation                                                                                                                               | Autoclave | Supplements |
|-----------------|---------------------------------------------------------|------------------------------|--------------------|-------------------------------------|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|-------------|
| <b>BS08408</b>  | COMPASS® cc Agar - Selective supplement                 | 10 vials<br>pack qsf 500 mL  | 3 years / 30 days* | 365                                 | 2-8°C                      | 30 at 2-8°C                                                                                                                                             | 0         | -           |
| <b>BS08508</b>  | COMPASS® Bacillus cereus - Selective supplement         | 500 vials<br>pack qsf 500 mL | 2 years            | 240                                 | 2-8°C                      | 30 at 2-8°C                                                                                                                                             | 0         | -           |
| <b>BS08908</b>  | RHAPSODY supplement                                     | 500 vials<br>pack qsf 500 mL | 3 years            | 365                                 | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS09008</b>  | EASY Staph® supplement                                  | vials qsf 100                | 3 years            | 365                                 | 2-8°C                      | 8 days at 2-8 °C.Reheat to between 25 and 37 °C before use. The white precipitate observed at 2-8 °C will disappear when heated to between 25 and 37°C. | -         | -           |
| <b>BS09208</b>  | D-cycloserine Liquid supplement                         | 10 vials qsf 90 mL           | 1 an               | 90                                  | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS09308</b>  | IRIS <i>Salmonella</i> ® supplement qsf 90 mL - Tablets | 120 tablets                  | 42 months          | 180                                 | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS09408</b>  | D-cycloserine Liquid supplement                         | 1 vial qsf 50 mL             | 1 an               | 90                                  | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS09508</b>  | CSD® supplement                                         | 10 × 100 mL in vials         | 1 year             | 90                                  | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS09708</b>  | EASY <i>Staph</i> ® Supplement                          | 1 vial qsf 500 mL            | 3 years            | 365                                 | 2-8°C                      | 8 days at 2-8 °C.Reheat to between 25 and 37 °C before use. The white precipitate observed at 2-8 °C will disappear when heated to between 25 and 37°C. | 0         | -           |
| <b>BS09808</b>  | COMPASS® Bacillus Plus Agar - Complete supplement       | 10 x 50 mL in vials          | 365 days           | 90                                  | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS09908</b>  | CSD® supplement qsf 10 g                                | 120 tablets                  | 24 mois            | 180                                 | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS10008</b>  | CSD® supplement qsf 25 g                                | 120 tablets                  | 24 mois            | 180                                 | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS10208</b>  | IRIS <i>Salmonella</i> ® concentrated supplement        | 10 x 90 mL in vials          | 18 months          | 90                                  | 2-8°C                      | -                                                                                                                                                       | -         | -           |
| <b>BS10408</b>  | PERFECTO Duo Agar liquid selective supplement           | al of 50 ml qs               | Not created yet    | pas encore créée                    | 2-8°C, shielded from light | -                                                                                                                                                       | -         | -           |

| Designation ENG | Packaging                                            | Packaging                                                                                                                                           | Shelf life           | Minimal<br>Remaining<br>shelf life<br>(days) | Storing T° | Storage after preparation              | Autoclave | Supplements |
|-----------------|------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|----------------------------------------------|------------|----------------------------------------|-----------|-------------|
| <b>BS10608</b>  | Oxford selective supplement                          | 10 vials<br>pack qsf<br>500 mL                                                                                                                      | 3 years / 7<br>days* | 365                                          | 2-8°C      | 7 days at 2-8°C shielded from<br>light |           |             |
| <b>BT00108</b>  | Microplate MUG/EC                                    | 25 opaque<br>microplates                                                                                                                            | 730 days             | 90                                           | 2-8°C      | -                                      | -         | -           |
| <b>BT00308</b>  | Microplate MUD/SF                                    | 25 opaque<br>microplates                                                                                                                            | 730 days             | 90                                           | 2-8°C      | -                                      | -         | -           |
| <b>BT00508</b>  | BAIRD-PARKER RPF Agar                                | 6 × 90 mL<br>vials of<br>base ready-<br>to-melt<br>+ 6 vials of<br>RPF<br>supplement<br>qsf 100 mL                                                  | 365 days             | 90                                           | 2-8°C      | -                                      | -         | -           |
| <b>BT00808</b>  | COMPASS® <i>Listeria</i> Agar                        | 6 × 200 mL<br>base<br>medium in<br>vials<br>+ 6 vials of<br>selective<br>suppl. qsf<br>200 mL<br>+ 6 vials of<br>enrichment<br>suppl. qsf<br>200 mL | 365 days             | 90                                           | 2-8°C      | -                                      | -         | -           |
| <b>BT01008</b>  | BAIRD-PARKER RPF Agar                                | 6 × 180 mL<br>vials of<br>base ready-<br>to-melt<br>+ 6 vials of<br>RPF<br>supplement<br>qsf 200 mL                                                 | 365 days             | 90                                           | 2-8°C      | -                                      | -         | -           |
| <b>BT01108</b>  | CONFIRM' <i>Salmonella</i> - Latex agglutination kit | 50 tests                                                                                                                                            | 18 months            | 90                                           | 2-8°C      | -                                      | -         | -           |



| Designation ENG   | Packaging                                   | Packaging                                                                                    | Shelf life | Minimal<br>Remaining<br>shelf life<br>(days) | Storing T° | Storage after preparation | Autoclave | Supplements |
|-------------------|---------------------------------------------|----------------------------------------------------------------------------------------------|------------|----------------------------------------------|------------|---------------------------|-----------|-------------|
| <b>BT01208</b>    | EASY <i>Staph</i> ® Agar                    | 6 × 190 mL<br>vials of<br>base ready-<br>to-melt<br>+ 6 vials of<br>supplement<br>qsf 200 mL | 365 days   | 90                                           | 2-8°C      | -                         | -         | -           |
| <b>BT01308</b>    | EASY <i>Staph</i> ® Agar                    | 6 × 90 mL<br>vials of<br>base ready-<br>to-melt<br>+ 6 vials of<br>supplement<br>qsf 100 mL  | 365 days   | 90                                           | 2-8°C      | -                         | -         | -           |
| <b>BT01408</b>    | <i>Listeria</i> - ID SYSTEM                 | 20 tests                                                                                     | 365 days   | 90                                           | 2-8°C      | -                         | -         | -           |
| <b>BT01508</b>    | <i>Legionella</i> - Latex agglutination Kit | 50 tests                                                                                     | 365 days   | 90                                           | 2-8°C      | -                         | -         | -           |
| <b>BT01708</b>    | COMFIRM'C. perfringens                      | qsf 50 tests<br>6 x R1 qsp,<br>2 x R2, 1 x<br>R3 +<br>barretes                               | -          | -                                            | 2-8°C      | -                         | -         | -           |
| <b>BT01708_R1</b> | COMFIRM'C. perfringens                      | qsf 50 tests<br>6 x R1 qsp,<br>2 x R2, 1 x<br>R3 +<br>barretes                               | -          | -                                            | 2-8°C      | -                         | -         | -           |
| <b>BT01708_R2</b> | COMFIRM'C. perfringens                      | qsf 50 tests<br>6 x R1 qsp,<br>2 x R2, 1 x<br>R3 +<br>barretes                               | -          | -                                            | 2-8°C      | -                         | -         | -           |
| <b>BT01708_R3</b> | COMFIRM'C. perfringens                      | qsf 50 tests<br>6 x R1 qsp,<br>2 x R2, 1 x<br>R3 +<br>barretes                               | -          | -                                            | 2-8°C      | -                         | -         | -           |